



## BSU810

### SUPPORTO DA PARETE OB WB

- Wall bracket anchorage for OB2012, OB3000 and OB1000 medical suction units
- Anchorage system which integrates the recharging system
- Anchorage system tested up to 10g capacity as obligatory as per EN 1789:1999 directive
- Extremely compact and functional wall bracket, not protruding out of rescue vehicle wall, with fast release button. Complete with cable for the connection to the ambulance electric system
- 12V power supply for OB2012 and OB1000; 11÷30V power supply for OB3000
- Weight: about 0,8 kg
- Dimensions: about 240 x 255 x 5 mm

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.