



BSU212

OB 1000 FM MEDICAL SUCTION UNIT

- Medical Suction Unit OB 1000 powerful, easy and compact
- Max. suction flow 25 liter/minute
- Max suction rate 850 mbar (80 kPa)
- Autonomy on continuous functioning >45 Minutes.
- New canister OB-J features one canister and one adapter that allows the use of the jar with the disposable liner produced by the company Serres
- Cable for the connection to the vehicle lighter (12Vcc)
- Can be provided with a 2-poles male connector main power supply, (power supply voltage rate between 100 e 240 Vac 50/60 Hz)
- Dimensions: about 335 x 290(h) x 105 mm
- Weight: about 3,3 kg
- Manufactured in compliance with all current norms

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.