



BSU368EU

OB MINI FM MEDICAL SUCTION UNIT BLU EU

- NEW Medical Suction Unit OB MINI
- Max. vacuum rate 850 mbar $\pm 10\%$
- Max. flow rate 25 liter/minute
- Lithium polymer battery (solid-non organic) allowing the autonomy of more than 70 minutes on free cycle
- New 1l canister OB-J for secretion collection made in polypropylene for medical use, the jar features one canister and one adapter that allows the use of the jar with the disposable liner produced by the company Serres (not included)
- Battery status indicator with no need to turn the device on
- Supplied with recharging cable (11÷25 Vcc)
- OB MINI features the inner microprocessor for the production data record and functioning tracing
- Safety inspection reminder indicates that the suction unit has to be sent for the scheduled maintenance
- Blue transport bag
- Supplied with universal power supply (from 90 to 240 Vac) EU, UK or JP/USA
- Compliant to all current norms, 93/42/CEE Directive and main reference norm ISO 10079-1
- Dimensions: 290x240x95 mm
- Weight: 1,8 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.